



Pride® Lift Chairs
Experience the Luxury

Pride®'s Classic Collection delivers superb value and excellent comfort with features like an innovative arm design and overstuffed fiber-filled backs. In addition, Classic Collection Lift Chairs feature patented, quiet and smooth operation and stylish designs to add to any home décor.

Classic COLLECTION

CLASSIC COLLECTION MODELS



CL-10	CL-15 ²	CL-20	CL-30 ²
2-Position Medium Button-back Design 5'2" - 5'10"	3-Position Medium Button-back Design 5'2" - 5'10"	2-Position Medium Sewn Pillow-back Design 5'4" - 6'0"	3-Position Medium Chaise Lounger Split-back Design 5'4" - 6'0"

SPECIFICATIONS

	CL-10	CL-15	CL-20	CL-30
Seat-To-Floor:	20"	20"	20"	20"
Seat Depth:	20"	20"	20"	20"
Seat Width:	20"	20"	20"	22"
Top of Back to Seat:	25"	27"	24.5"	26"
Distance From Wall:	12"	18"	11"	18"
Overall Chair Width:	33.5"	33.5"	32.5"	34"
Chair Weight:	110 lbs.	119 lbs.	114 lbs.	118 lbs.

BENEFITS

- Engineered furniture grade laminate/ hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head & arm covers
- Fire retardant (California specifications; bulletin 117)
- Standard, large dual pockets
- Removable back - no tools required
- 375 lbs. weight capacity on CL-20 & CL-30; 325 lbs. weight capacity on CL-10 & CL-15
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 12' cord
- Lifetime warranty on steel lift frame
- 2-year limited warranty

FABRIC SELECTION

STANDARD FABRICS

MARINE WINE MOSS WHEAT TAUPE

OPTIONAL FABRICS

RASPBERRY HERB ALMOND TRUFFLE BRASS ROSE TWIST OATMEAL SLATE

PADDED-SUEDE

BEACH OLIVE BLUE EXPRESSO

NEW LEXIS VINYL

BLACK MUSHROOM BURGUNDY NEW CHESTNUT ROYAL BLUE

PLUSH VELVET COLLECTION

BERRY CHOCOLATE COFFEE

WINDSOR MICRO-SUEDE COLLECTION

CAMEL MAPLE

NEW CRYPTON® SUPER FABRICS

SMART SUEDE**

OYSTER PICKLE CLOVER MALT GINGER GUMBALL SOY TERRA COTTA CAPPUCCINO

SMART FABRICS**

VANILLA BEAN RED PEPPER CELERY CREAM

CRYPTON®
Super Fabrics


OPTIONS*

- Heat & Massage (CL-20 & CL-30 only)
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers
- Leather/vinyl wipes
- Footrest extension available on CL-15 & CL-30

Additional lead-time may apply for some options.

* Additional charges apply for all options.

1. Batteries not included.
2. Available with optional footrest extension.



Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.
**5 year warranty.