



THE CITY LIFE.

With sleek good looks and high caliber performance, the C300 Corpus offers the perfect blend of stability and agility. Features include: adjustable suspension; front and rear battery access; and a seat-to-floor height of 17.5 inches with 45 degrees of tilt. No other product in its class can even compete. Small, nimble and quick, the C300 Corpus is designed for navigating busy streets and getting in and out of tight spaces. Plus, the Corpus is the world's most popular powered seating system with an ergonomic back featuring built-in, adjustable lumbar and trunk supports. The modular design of the Corpus also allows for growth and added seat functions. An ideal pairing, the C300 Corpus is ready to take you in any direction you choose. Lead the way.

C300
CORPUS



POWER FUNCTIONS



Power adjustable tilt



R-Net electronics with ICS



Power adjustable recline










Power seat height adjustment



Power adjustable legrests

SHROUD COLORS

-  Crystal Blue
-  Vintage Wine
-  Metallic Green
-  Platinum Silver
-  Black Onyx
-  Chili Pepper Red
-  Popstar Pink
-  Titan Blue

SPECIFICATIONS

Maximum Speed	5 MPH	Gradient Capacity	10° Up 6° Dn
Driving Range	15-19 Miles	Power Seat Elevator	8"
Max User Weight	265 lbs	Seat Widths	17" - 19" (Adjustable)
Ground Clearance	3"	Seat Depths	16" - 22" (Adjustable)
Battery Capacity	Group 34	Backrest Heights	16" - 22"
Width	24.25" - 28.5"	Seat To Floor Height	17.25" - 21.5"
Length	43"	Electronics (standard)	R-Net
Recommended Battery	Group 34	Alternative Drive Controls	Yes w/R-Net



THE POWER OF MOBILITY

Permobil Inc.
 300 Duke Drive
 Lebanon, Tennessee 37090
 800.736.0925
 www.permobil.com