



WRANGLER®

The Wrangler® PMV is the ultimate personal mobility vehicle, providing optimal power, superior performance, unmatched reliability and superb style.

- Dual motors for maximum traction and power
- 13" all-terrain tires for enhanced performance on rugged terrain
- Deluxe contour high-back seat with headrest and sliders
- Standard deluxe contoured high-back seat
- 400 lbs. weight capacity and up to 10 mph top speed



PURSUIT®

The Pursuit® defines comfort and safety with a reclining high-back seat with headrest and sliders, a full lighting package with directional signals and lower headlight, and a standard front basket and rearview mirror.

- Deluxe reclining high-back seat with headrest and sliders
- Full lighting package including directional signals and angle-adjustable, lower headlight
- Front basket and rearview mirror standard
- Large 13" low-profile solid tires for excellent outdoor performance
- 400 lbs. weight capacity and up to 7.25 mph top speed



PURSUIT® XL

The Pursuit® XL is an exciting blend of power and precision that features full suspension, low-profile tires and a hydraulic sealed brake system, making it the ultimate outdoor scooter.

- Front and rear suspension
- Large 13" low-profile solid tires for excellent outdoor performance
- Features 120 amp microprocessor based controller for powerful drivetrain
- Full lighting package with directional signals and angle-adjustable, lower headlight
- 400 lbs. weight capacity and up to 8.25 mph top speed



www.pridemobility.com

800-800-8586 (US) • 888-570-1113 (Canada)

Pride® Scooters

The Ultimate in Style & Performance





GO-GO® ULTRA X

The Go-Go® Ultra X is loaded with design features like one-hand feather-touch disassembly and a convenient drop-in battery box that makes transport and travel worry-free.

- Auto-latching lockup mechanism permits simple frame separation with only one hand
- All-new frame design easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Available in 3 or 4-wheel format
- 260 lbs. weight capacity and up to 4 mph top speed



GO-GO ELITE TRAVELLER®

The Go-Go Elite Traveller®'s compact design allows it to easily maneuver in tight spaces while providing stable outdoor performance and one-hand feather-touch disassembly.

- Auto-latching lockup mechanism permits frame separation with only one hand
- All new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Includes two sets of easily changeable color panels – Red and Blue
- Available in 3 or 4-wheel format
- 275 lbs. weight capacity and up to 4.25 mph top speed



GO-GO ELITE TRAVELLER® PLUS/PLUS HD

The Go-Go Elite Traveller® Plus offers increased length, width, and weight capacity and a wraparound delta tiller to allow larger individuals and those with limited dexterity use of a travel mobility product.

- Delta tiller with wraparound handles
- Standard 18AH battery pack for enhanced range
- Includes two sets of easily changeable color panels – Red and Blue
- Available in 3 or 4-wheel format
- 300 lbs. weight capacity (350 lbs. weight capacity Plus HD version) and up to 4.5 mph top speed

VICTORY® 9

The Victory 9® delivers the high performance operation and feature-rich design of the Victory family of scooters.

- Full complement of advanced standard features including a wraparound delta tiller
- 9-inch wheels and compact overall size for enhanced maneuverability
- Front running light
- Available in 3 or 4-wheel format
- 300 lbs. weight capacity and up to 5.25 mph top speed



VICTORY® 9 PS

The Victory 9® PS 3-wheel delivers a long list of exciting features including a power elevating seat, low-profile tires, feather-touch disassembly, backlit battery gauge and wraparound delta tiller.

- Standard power elevating seat
- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Delta tiller with wraparound handles
- Long-lasting, high intensity LED headlight offers optimal pathway illumination
- 300 lbs. weight capacity and up to 5.5 mph top speed



VICTORY® 10

The Victory® 10 combines superb comfort and solid performance with luxury scooter touches at a great low price. It features a standard front basket and headlight.

- Full complement of advanced standard features
- Pride's exclusive low-profile, black, non-scuffing wheels
- Delta tiller with wraparound handles
- Available in 3 or 4-wheel format
- 400 lbs. weight capacity and up to 5.25 mph top speed



CELEBRITY® X

Innovative features, including a delta tiller and bright long lasting LED lights, make the Celebrity® X the best overall scooter value in its class.

- Feature-rich luxury product at an attractive price
- Delta tiller with wraparound handles
- Auto-latching lockup for easy assembly/disassembly
- Available in 3 or 4-wheel format
- 350 lbs. weight capacity and up to 6 mph top speed

MAXIMA

The new, restyled Maxima delivers luxury and heavy-duty performance at an exceptional value. Engineered with the larger individual in mind, the Maxima features a weight capacity of 500 lbs.

- Standard conveniences and safety features like a front basket and full lighting package
- Ultra heavy-duty drivetrain with 110-amp controller
- Ergonomically enhanced 22" wide x 18.5" deep seating
- Available in 3 or 4-wheel formats
- 500 lbs. weight capacity and up to 5.25 mph top speed

